

Characterization of Respirable Output from the Vyntus™ APS dosimeter used with the Cirrus2 nebulizer

Author: Margriet Kaldenbach

Jaeger Medical

Introduction

With the introduction of the new Cirrus™2 nebulizer, the Vyntus™ APS system was characterized to assess the respirable output of this configuration. This white paper summarizes the methods and results, providing health care professionals with practical data to support consistent bronchial challenge testing (BCT) in daily practice.

Consistent and well-controlled aerosol delivery is important for BCT. Various aerosol delivery systems are available, including jet nebulizers, breath-actuated nebulizers, and dosimeter-based devices. The Vyntus™ APS is a dosimeter device, delivering short, precisely timed aerosol actuations (0.2–1.0 seconds) during inhalation, ensuring consistent dosing for each breath. Whenever components of a dosimeter system are updated, their performance should be re-evaluated to maintain consistent delivered doses and reproducible BCT results over time.^{1,2}

The 2017 ERS/ATS technical standards for methacholine challenge testing¹ recommend that devices intended for BCT be characterized to confirm respirable output. This includes measuring inhaled mass (IM) via filter-based aerosol collection and determining the respirable fraction (RF; aerosol droplets $\leq 5 \mu\text{m}$) under simulated breathing conditions.^{1,2} For a dosimeter such as the Vyntus™ APS, characterization is performed using actuation times of 0.2–1.0 seconds, reflecting typical clinical practice.¹ The product of IM and RF defines the respirable output, a key parameter for calculating delivered doses and setting up challenge protocols.

Methods characterization Vyntus™ APS with Cirrus2 nebulizer

A series of tests was conducted to evaluate the respirable output in mg/min of the newly introduced Cirrus™2 nebulizer when paired with the Vyntus™ APS dosimeter. The procedures below describe the system setup and how IM and RF were determined.

Vyntus™ APS system setup

One Vyntus™ APS system was used for testing, including the APS Connection Block for Cirrus™2 (V-707506, Jaeger Medical GmbH) and the APS Tube (V-892080, single-use, Jaeger Medical GmbH), see


Figure 1. Five Cirrus™2 nebulizers (V-8922609, single-use, Intersurgical® GmbH, Germany) from three production lots were tested, each paired with its own dedicated APS Intersurgical® filter (V-892060, single-use, Intersurgical® GmbH, Germany). These pairings were maintained throughout all tests. Prior to testing, all five nebulizer/filter sets underwent a compressed air check, all checks were within the limit range of 0.09–0.170 L/s. Testing was conducted using SentrySuite™ software version 3.30.



Figure 1. Vyntus™ APS system set up

Test procedure Inhaled Mass

All five nebulizer/filter sets were evaluated under controlled measurement conditions (see **Figure 2**). Each nebulizer was filled with 3 mL of 0.0625% albuterol hemisulfate solution containing 0.5% fluorescein, allowing measurement via



ultraviolet–visible (UV-VIS) spectroscopy. Aerosol output was collected on inhalation filters (V-892020, AirLife™, single-use, Vyaire Medical, USA) at the patient outlet, with a breathing simulator (BRS 1100,



Copley Scientific, UK) connected to reproduce a 1:1 inhalation-to-exhalation ratio, tidal volume 750 mL, and 15 breaths per minute. The Vyntus™ APS was tested at three actuation durations (0.2, 0.4, 1.0 seconds) with 30, 15, and 10 actuations per setting, respectively. Each nebulizer was tested five times per setting, yielding 75 IM measurements in total. Inhaled mass ($\mu\text{g}/\text{sec}$) was determined by UV–VIS spectroscopy. Inhaled volume ($\mu\text{L}/\text{sec}$) was calculated by dividing recovered dose by total nebulization time, then converted to IM (mg/min) using the formulation density ($1.0001 \text{ mg}/\mu\text{L}$) and converting seconds to minutes.

Figure 2. System setup for IM testing

Test procedure Respirable fraction

RF was measured using the same nebulizer/filter sets and breathing simulator settings (see **Figure 3**). Each nebulizer was filled with 3 mL of 0.0625% albuterol hemisulfate solution (without fluorescein). Aerosol particle size

distribution was determined by laser diffraction (Malvern Spraytec) to calculate the RF. The same actuation durations (0.2, 0.4, 1.0 seconds) and actuation counts (30, 15, 10) were applied, with five tests per nebulizer per setting, resulting in 75 RF measurements.

Figure 3. System setup for RF testing
Results Characterization Cirrus™2 nebulizer when used with Vyntus™ APS Dosimeter

Results for both IM and RF are summarized in **Table 4**. For each actuation time (0.2, 0.4, and 1.0 seconds), the table presents mean values and the corresponding coefficient of variation (CoV in %) across all five nebulizers. It also shows the derived respirable output (mg/min), calculated as the product of inhaled mass and respirable fraction for each actuation time.

	Actuation time		
	0.2 sec	0.4 sec	1.0 sec
Inhaled mass in mg/min (CoV)	243 (9.7)	262 (10.4)	275 (8.5)
Respirable fraction in % (CoV)	74 (4.1)	75 (5.1)	80 (4.0)
Respirable output in mg/min (CoV)	180 (10.5)	197 (11.6)	220 (9.4)

Table 4. Results IM, IF and Respirable Output

The overall mean respirable output across the three actuation times was

approximately 200 mg/min. Considering the inherent variability of BCT testing, this rounded value was selected as a representative estimate for respirable output and can be applied to actuation times between 0.2 and 1.0 seconds.

Conclusion

Characterization of the Vyntus™ APS Dosimeter with the new Cirrus™2 nebulizer confirmed that it delivers consistent aerosol performance under controlled conditions. Both IM and RF showed stable results across actuation times of 0.2, 0.4, and 1.0 seconds, resulting in a mean respirable output of approximately 200 mg/min. This value provides a practical basis for developing standardized BCT protocols, allowing the delivered dose to be accurately calculated at each step. Considering the inherent variability in BCT testing^{1,3,4} the rounded mean of 200 mg/min can be applied across the full actuation time range of 0.2–1.0 seconds, ensuring consistent dose delivery. These reference data support clinical users in confidently applying the device with nebulizers for methacholine challenge testing, in line with current ERS/ATS 2017 recommendations.¹

References

1. A. L. Coates et al., "ERS technical standard on bronchial challenge testing: General considerations and performance of methacholine challenge tests," *Eur. Respir. J.*, vol. 49, no. 5, May 2017, doi: 10.1183/13993003.01526-2016.
2. A. Coates et al., "Characterizing nebulizer performance for methacholine challenge tests," *Am. J. Respir. Crit. Care Med.*, vol. 198, no. 8, pp. 988–990, 2018, doi: 10.1164/rccm.201805-0942ED.
3. R. O. Crapo et al., "American Thoracic Society guidelines for methacholine and exercise challenge testing -1999," *Am. J. Respir. Crit. Care Med.*, vol. 161, no. 1, pp. 309–329, 2000.
4. T. S. Hallstrand et al., "ERS technical standard on bronchial challenge testing: Pathophysiology and methodology of indirect airway challenge testing," *Eur. Respir. J.*, vol. 52, no. 5, pp. 1–18, 2018, doi: 10.1183/13993003.01033-2018.

This document is intended for health care professionals for educational purposes only.



 Jaeger Medical GmbH
Leibnizstrasse 7
97204 Hoechberg
Germany

Jaeger Medical America, Inc.
14050 Rebecca St, Bldg 4
Moreno Valley, CA 92553
USA

For International distribution

© 2025 Jaeger Medical GmbH. Jaeger and the Jaeger logo are trademarks or registered trademarks of Jaeger Medical GmbH or one of its affiliates. Jaeger Medical devices are class I & IIa according to Medical Devices Directive 93/42/EEC or Medical Device Regulation EU 2017/745 as indicated on each declaration of conformity. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labeling. JAE-INT-2500053 1.0