



Return-Shipment of Parts containing Lithium-Ion-Batteries

For a secure return shipment of Parts containing Lithium-Ion-Batteries, the following instructions must be followed carefully:

Packaging:

Whenever possible, please use the original packaging with carton and foam inlays. When that is not possible, please use a stable carton and enough packaging material to ensure that the appliance is secured.

Marking:

When using the original packaging, please check that the labels are in good condition and readable. If this is not possible, or you used your own carton, please ensure to use the necessary label (UN3481, Part-nr. V-919840) The minimum dimensions of the label are 120mm x 110mm.

Documentation:

- for Airfreight-Shipments, the statement "**Lithium-ION batteries not restricted as per PI 967 Part II** " has to be added on the Airwaybill (AWB)

What to do if the battery is damaged?

DO NOT SHIP DAMAGED BATTERIES! It is not permitted by law. Damaged batteries must be disposed of at your local recycling depot.

IMPORTANT NOTE: The sender is responsible for safely and correctly packing and labeling the package. The sender is also responsible for the documentation and declaration of the lithium-ion-battery according to the applicable dangerous goods laws.

If you have any questions, please contact Jaeger Medical at shipping.eu@jaegerrdx.com